

## Product Information Sheet

### Enabling Devices: Auditory Communicator



#### Enabling Devices Contact:

Telephone: 800-832-8697

Website: [www.enablingdevices.com](http://www.enablingdevices.com)

#### Support & Training:

- For general support and resources visit: <https://enablingdevices.com/resources>
- Enabling Devices Auditory Communicator User Guide: <https://enablingdevices.com/files/downloads/f2a11b79-6c0f-44a2-b5b9-917fa50ae056.pdf>

#### Auditory Communicator Features:

Display Type	Sequenced recorded messages with switch plate.
Dimensions	6½"L x 4"W x 2½"H
Weight	0.5 lbs.
Type of Voice Output	Recorded speech
Language System	The auditory communicator will scan through your pre-recorded messages at a lowered volume. When you reach the desired message, pushing the switch plate again plays the chosen message in an amplified volume. 16-level unit with 300 seconds of total recording time with 18 seconds per level.
Access Methods	In addition to a built in switch, there is a jack for an external switch (not included).
Battery	Requires 4 AA Batteries
Available Colors	Black, Blue