

# Product Information Sheet

## Enabling Devices: Tactile Symbol Communicator



### Enabling Devices Contact:

Telephone: 800-832-8697

Website: [www.enablingdevices.com](http://www.enablingdevices.com)

### Support & Training:

- For general support and resources visit: <https://enablingdevices.com/resources>
- Enabling Devices Communicator for the Blind User Guide:  
<http://enablingdevices.com/files/downloads/8859c4e6-aca3-4f51-b6e2-ecc01b51f323.pdf>
- Product video:  
[http://enablingdevices.com/catalog/products\\_for\\_the\\_visually\\_impaired/communicators\\_for\\_the\\_visually\\_impaired/tactile-symbol-communicator](http://enablingdevices.com/catalog/products_for_the_visually_impaired/communicators_for_the_visually_impaired/tactile-symbol-communicator)

### Tactile Symbol Communicator Features:

Display Type	6 separate 3" x 3" cells with 1" black separation between cells. Can be used with tactiles and/or small, lightweight objects.
Dimensions	14" L x 10 1/2" W x 2"H
Weight	3.5 lbs. with battery
Type of Voice Output	Recorded voice
Language System	Use with tactile symbols. Starter kit includes symbols for bathroom, something wrong, drink, I want, finished and eat. Stores 6 messages on 6 levels for a total of 36 messages. The vibration setting can be turned on or off. This feature may be useful for those who have a dual sensory impairment.
Battery	Requires 4 AA Batteries
Access Methods	Direct selection using light push pressure. Optional vibration feature.
Available Colors	Black