

# Product Information Sheet

## Prentke-Romich Company (PRC): Accent 1400



### PRC Contact:

(800) 262-1990

Hours M-F, 8:00 a.m. to 7:00 p.m. EST

Saturday-Sunday, 1:00 p.m. to 5:00 p.m. EST

[service@prentrom.com](mailto:service@prentrom.com)

### Support & Training:

- **PHX and Northern AZ:** Contact Jeremy Legaspi by e-mail at [Jeremy.Legaspi@prentrom.com](mailto:Jeremy.Legaspi@prentrom.com) or by phone: (800) 262 1984 ext. 452 or 602-717-5506 (cell).
- **Tucson and Southern AZ:** Contact LaConda Walker by e-mail at [LaConda.Ross@prentrom.com](mailto:LaConda.Ross@prentrom.com), or by phone: 1-800-484-8008 or 520-870-4648.
  - <https://www.prentrom.com/support> includes downloads, live chat, and upgrade information. PASS is also available for all Accent products and can be downloaded at no cost. Therapy materials/pictures can be generated to engineer environments, modify books, and create learning materials using this program.
  - PRC Accent 1400 Quick Reference PDF:  
<https://www.prentrom.com/support/accent/download-accent-1400-getting-started-2> or  
<http://file.prentrom.com/278/17372v1.3-Accent%201400%20Getting%20Started.pdf>
  - Product training:  
[https://www.prentrom.com/training/classes?type=master the technology&method=online&page=1](https://www.prentrom.com/training/classes?type=master%20the%20technology&method=online&page=1)
  - Implementation Class :  
[https://www.prentrom.com/training/classes?type=making aac work&method=online&page=1](https://www.prentrom.com/training/classes?type=making%20aac%20work&method=online&page=1)

- Live and On Demand training: <https://www.prentrom.com/training/>
- You Tube Channel for programming Accent Devices: <https://www.youtube.com/user/PRCaccess/playlists>
- Learn about Minspeak language system <http://minspeak.com/>. The “what,” “who,” and “why” of Minspeak.
- PRC Language Lab at <http://www.aaclanguagelab.com/resources>. Access to this site, via paid subscription (\$84), provides many different materials related to language development for AAC users.
- Realize Language: <https://realizelanguage.com/info/faqs> Subscription tool (\$84 annually) to analyze device user’s messages.

### Accent 1400 Features:

The cost to unlock this device with PRC is \$50

Display Type	14.1“ Dynamic display touchscreen. Manufactured with chemically resistant glass. Screen Resolution 1280 x 800
Dimensions	14.1" w × 8.9" h × 1.8" d
Weight	5.9 lbs. with battery
Type of Voice Output	Synthesized Speech with Digitized Speech Capacity. Voices include: Nuance, Microsoft, Ivona American English and Ivona American Spanish
Language System	The Accent 1400 offers a choice of vocabulary options featuring Unity® with LAMP/Words for Life™, CoreScanner™, and WordCore (allowing users to generate spontaneous and specific messages). Other available vocabulary options include Essence™ for literate adults and WordPower™
Operating System	Windows® operating system with Bluetooth™ and WiFi powered by a Multi-Core Bay Trail Intel™ processor. The device can be unlocked to access Windows applications and integrated to allow Wi-Fi™ connectivity.
Battery Run Time/Charger	18 to 20 hours depending on individual usage and power settings; however, while using NuEye battery life is 10 to 12 hours.
Access Methods	Multiple access options include direct selection, 1 and 2 switch scanning, mouse emulation, joystick, integrated head mouse with NuPoint, and integrated eye gaze with NuEye.
Other	3 - 2.0 USB ports and 1 - 3.0 USB port, 1 - HDMI port, carrying handle, front and rear cameras, and Integrated Bluetooth® connectivity
Available Colors	Black, Blue, Green, Red, Yellow