

Product Information Sheet

Prentke-Romich Company (PRC): Accent 800



PRC Contact:

(800) 262-1990

Hours M-F, 8:00 a.m. to 7:00 p.m. EST

Saturday-Sunday, 1:00 p.m. to 5:00 p.m. EST

service@prentrom.com

Support & Training:

- **PHX and Northern AZ:** Contact Jeremy Legaspi by e-mail at Jeremy.Legaspi@prentrom.com or by phone: (800) 262 1984 ext. 452 or 602-717-5506 (cell).
- **Tucson and Southern AZ:** Contact LaConda Walker by e-mail at LaConda.Ross@prentrom.com, or by phone: 1-800-484-8008 or 520-870-4648.
 - <https://www.prentrom.com/support> includes downloads, live chat, and upgrade information. PASS is also available for all Accent products and can be downloaded at no cost. Therapy materials/pictures can be generated to engineer environments, modify books, and create learning materials using this program.
 - PRC Accent 800 Quick Reference PDF: <https://www.prentrom.com/downloads/device-and-accessory-manuals?id=251>
 - Implementation Class :
https://www.prentrom.com/training/classes?type=making_aac_work&method=online&page=1
 - Live and On Demand training: <https://www.prentrom.com/training/>

- You Tube Channel for programming Accent Devices:
<https://www.youtube.com/user/PRCaccess/playlists>
- Learn about Minspeak language system <http://minspeak.com/>. The “what,” “who,” and “why” of Minspeak.
- PRC Language Lab at <http://www.aaclanguagelab.com/resources>. Access to this site, via paid subscription (\$84), provides many different materials related to language development for AAC users.
- Realize Language: <https://realizelanguage.com/info/faqs> Subscription tool (\$84 annually) to analyze device user’s messages.

Accent 800 Features:

The cost to unlock this device with PRC is \$50

Display Type	8” Dynamic touchscreen display. Screen Resolution 1280 x 800
Dimensions	9" w × 5.8" h × 1.4" d
Weight	2.4 lbs.
Type of Voice Output	Synthesized Speech with Digitized Speech Capacity. Voices include: Nuance, Microsoft, Ivona American English and Ivona American Spanish.
Language System	The Accent 800 offers a choice of vocabulary options featuring Unity® with LAMP/Words for Life™, CoreScanner™, and WordCore (allowing users to generate spontaneous and specific messages). Other available vocabulary options include Essence™ for literate adults and WordPower™
Operating System	Microsoft Windows®; The device can be unlocked to access Windows applications and integrated to allow Wi-Fi™ connectivity.
Battery Run Time/Charger	11 to 13 hours depending on individual usage and power settings
Access Methods	Multiple access options include direct selection, 1 and 2 switch scanning (integrated ports), mouse, and joystick
Other	Includes a carrying case for portability, front and rear cameras, Integrated Bluetooth® connectivity, and IR controls. 1 - 3.0 USB port and 1 - 2.0 USB port, 1 mini USB port
Available Colors	Black, Blue, Red, Green, Yellow, Purple and Pink