

Product Information Sheet Saltillo Inc: TouchChat Express



Saltillo Contact:

Telephone: 1-800-382-8622

Email: info@saltillo.com

Support & Training:

- **PHX and Northern AZ:** Contact local Saltillo representative Sheri Predebon by e-mail at sheri@saltilloconsultant.com or by phone at 602-295-9987.
- **Tucson and Southern AZ:** Contact local Saltillo representative Jane Anderson by e-mail at jane@saltilloconsultant.com or by phone at 520-329-0624.
- Live and recorded webinars and video training at <https://saltillo.com/training>.
- For Live Chat, support articles and downloads, visit: <https://saltillo.com/support>. Live Chat is available Monday-Friday 9:00 am – 5:00 pm ET.
- To speak with a technician contact: 1-800-382-8622 or email: service@saltillo.com
- TouchChat Express Reference Guide: <https://touchchatapp.com/support/articles/faq/quick-reference-guide>

TouchChat Express Features:

Display Type	Touchscreen with retina HD display iPad Air: 9.7" iPad Mini: 7.9"
Dimensions	iPad Air: 6.88" x 9.63" x 1.68" iPad Mini: 5.50" x 8.13" x 1.68"
Weight	iPad Air: 1.90 lbs. iPad Mini: 1.35 lbs.
Type of Voice Output	Words, phrases and messages are spoken with a built-in Acapela voice synthesizer or by playing a recorded message. Various English and Spanish synthesized voices are available.
Language System	TouchChat Express includes 4 vocabulary packages: VocabPC, Multichat 15, Spelling (with 4 cell word prediction) & Primary. It includes both English and Spanish vocabularies, including WordPower and iEssence. Gateway and Pixon are available for purchase. All page sets can be modified and customized.
Operating System	iOS
Other	Front and rear camera
Battery Run Time/Charger	Up to 10 hours. Rechargeable.

If the end user of this device has a ventricular shunt, please note that some devices such as the iPad 2nd and 3rd generation versions contain magnets inside them. There have been studies done indicating that magnets may alter ventricular shunt settings if placed in very close contact (within 10 cm) to the programmable valve of the shunt. You can minimize the risk of this happening by keeping toy magnets, electronic devices, and any other devices containing magnets at a safe distance away from the head. When regularly using any of these devices, keep them from coming in close contact (within 10 cm) to the valve. The original version of the iPad does not contain magnets, so this caution does not apply to first generation iPads.