

# Product Information Sheet

## Tobii Dynavox: T15



### Tobii Dynavox Contact:

Telephone: 1-800-344-1778

Website: <http://www.tobiidynavox.com/>

Email: [support.na@tobiidynavox.com](mailto:support.na@tobiidynavox.com)

### Support & Training:

- Contact local Tobii DynaVox representative John Weber by e-mail at [john.weber@tobiati.com](mailto:john.weber@tobiati.com) or by phone at 480-220-8268.
- For general support visit: <http://www.tobiidynavox.com/support-training/>
- To helpful training videos visit: <https://www.youtube.com/user/DynaVoxVideos/videos> and <https://mydynavox.com/Resources/Videos>
- For support and downloads: <http://www.tobiidynavox.com/support/t-series/>
- Tobii Dynavox T15 User Manual: <https://fccid.io/document.php?id=2531515>
- Ideas for everyday use and therapy: <https://mydynavox.com/Resources/IdeasForTherapy>

### Tobii Dynavox T15 Features:

The cost to unlock this device with Tobii Dynavox is \$40.

Display Type	Ruggedized, backlit, high-definition display with touch screen
Dimensions	Screen size: 15.4" Device dimensions: 15.06" x 10.38" x 1.63 "
Weight	3 lbs. 13oz.
Operating System	Android
Type of Voice Output	Synthesized Acapela and IVONA voices for English. IVONA speech synthesizer includes adult, child and teen voices.
Language System	Compass™ is based on InterAACT, a communication and language framework. Additional options include Gateway, PODD, and WordPower page sets.
Access Methods	Touchscreen, keyboard, mouse/trackball, integrated single and dual switch scanning with audio feedback, head tracking, and keyguards are available.
Battery	Lithium Ion, more than 10 hrs. of typical use
Other	Device includes a rear and front facing camera, microUSB port, full size USB port, microHDMI port, headphone/microphone jack, and IR transceiver