

# Product Information Sheet

## Tobii Dynavox: i12



### Tobii Dynavox Contact:

Telephone: 1-800-344-1778

Website: <http://www.tobiidynavox.com/>

Email: [support.na@tobiidynavox.com](mailto:support.na@tobiidynavox.com)

Support hours: Weekdays, 9 am – 9 pm (US Eastern Standard Time)

### Support & Training:

- Contact local Tobii DynaVox representative John Weber by e-mail at [john.weber@tobiati.com](mailto:john.weber@tobiati.com) or by phone at 480-220-8268.
- For general support visit <http://www.tobiidynavox.com/support-training/>
- To view YouTube training videos visit: <https://www.youtube.com/channel/UCSFhjTdkb4rt4GilPmDypQA>
- i-series support page: <http://www.tobiidynavox.com/support/i-series-plus/>
- Tobii Dynavox i-series user manual: [http://www.toby-churchill.com/assets/downloads/Tobii\\_I-Series\\_User%E2%80%99s\\_manual\\_v\\_1\\_3\\_-\\_en-US\\_WEBHQ.pdf](http://www.toby-churchill.com/assets/downloads/Tobii_I-Series_User%E2%80%99s_manual_v_1_3_-_en-US_WEBHQ.pdf)

### Tobii Dynavox i12 Features:

The cost to unlock this device with Tobii Dynavox is \$40.

Display Type	LED backlit, dynamic display touchscreen. Gorilla® Glass, 10 point multi touch
Dimensions	Screen size: 12.1” Device dimensions: 12.1” x 10.8” x 4.1”
Weight	6.2 lbs
Operating System	Microsoft Windows
Type of Voice Output	Synthesized and recorded speech. Two built-in, closed box speakers. High quality sound at loud levels.
Language System	Includes Tobii Communicator software. Pre-made pages include text, single word, and phrase-based language. Custom page sets are often necessary and can be programmed. Compass, Gateway, and PODD vocabularies are available for additional purchase.
Access Methods	Built in touch screen, integrated single and dual switch scanning with auditory feedback, mouse, keyboard, and an optional IS20 Eye Tracker Module.
Battery Run Time	8 hours with eye tracker, 9.0 to 10.5 hours without eye tracker.
Other	Forward and rear facing cameras, microSD, Bluetooth