

Product Information Sheet

Tobii Dynavox: i12



Tobii Dynavox Contact:

Telephone: 1-800-344-1778

Website: <http://www.tobiidynavox.com/>

Email: support.na@tobiidynavox.com

Support hours: Weekdays, 9 am – 9 pm (US Eastern Standard Time)

Support & Training:

- Contact local Tobii DynaVox representative John Weber by e-mail at john.weber@tobiati.com or by phone at 480-220-8268.
- For general support visit <http://www.tobiidynavox.com/support-training/>
- To view YouTube training videos visit: <https://www.youtube.com/channel/UCSFhjTdkb4rt4GilPmDypQA>
- i-series support page: <http://www.tobiidynavox.com/support/i-series-plus/>
- Tobii Dynavox i-series user manual: http://www.toby-churchill.com/assets/downloads/Tobii_I-Series_User%E2%80%99s_manual_v_1_3_-_en-US_WEBHQ.pdf

Tobii Dynavox i12 Features:

The cost to unlock this device with Tobii Dynavox is \$40.

| | |
|----------------------|---|
| Display Type | LED backlit, dynamic display touchscreen. Gorilla® Glass, 10 point multi touch |
| Dimensions | Screen size: 12.1" Device dimensions: 12.1" x 10.8" x 4.1" |
| Weight | 6.2 lbs |
| Operating System | Microsoft Windows |
| Type of Voice Output | Synthesized and recorded speech. Two built-in, closed box speakers. High quality sound at loud levels. |
| Language System | Includes Tobii Communicator software. Pre-made pages include text, single word, and phrase-based language. Custom page sets are often necessary and can be programmed. Compass, Gateway, and PODD vocabularies are available for additional purchase. |
| Access Methods | Built in touch screen, integrated single and dual switch scanning with auditory feedback, mouse, keyboard, and an optional IS20 Eye Tracker Module. |
| Battery Run Time | 8 hours with eye tracker, 9.0 to 10.5 hours without eye tracker. |
| Other | Forward and rear facing cameras, microSD, Bluetooth |