

## Product Information Sheet

### Prentke-Romich Company (PRC): PRiO and PRiO-mini



#### PRC Contact:

(800) 262-1990

Hours M-F, 8:00 a.m. to 7:00 p.m. EST

Saturday-Sunday, 1:00 p.m. to 5:00 p.m. EST

[service@prentrom.com](mailto:service@prentrom.com)

#### Support & Training:

- **PHX and Northern AZ:** Contact Jeremy Legaspi by e-mail at [Jeremy.Legaspi@prentrom.com](mailto:Jeremy.Legaspi@prentrom.com), or by phone: (800) 262 1984 ext. 452 or 602-717-5506 (cell).
- **Tucson and Southern AZ:** Contact LaConda Walker by e-mail at [LaConda.Ross@prentrom.com](mailto:LaConda.Ross@prentrom.com), or by phone: 1-800-484-8008 or 520-870-4648.
  - <https://www.prentrom.com/support> includes downloads, live chat, and upgrade information. PASS is also available for Accent products and can be downloaded at no cost. Therapy materials/pictures can be generated to engineer environments, modify books, and create learning materials using this program.
  - PRiO-LAMP Words for Life Quick Reference PDF: <https://www.prentrom.com/support/attachments/280>
  - Installing and Setting up LAMP Words for Life on an iOS Device PDF: <https://www.prentrom.com/support/attachments/280>
  - Setting Up Your iPad for AAC: <https://aacapps.com/lamp/support/video/setting-up-your-ipad-for-aac>

- Live and On Demand training: <https://www.prentrom.com/training/>
- Realize Language: <https://realizelanguage.com/info/faqs> Subscription tool (\$84 annually) to analyze device user's messages.

## PriO Features:

Display Type	iPad Mini- 7.9 inch LED-backlit Multi-Touch display with IPS technology, 2048 × 1536 resolution with 326 pixels per inch iPad- 9.7 inch LED-backlit Multi-Touch display with IPS technology, 2048 × 1536 resolution with 3.1 million megapixels
Dimensions	iPad Mini- 7.87" (h) x 5.3" (w) x 0.29" (d) iPad- 9.4" (h) x 6.6" (w) x 0.24" (d)
Weight	iPad Mini- 0.73 lbs. iPad Air 2- 0.96 lbs.
Type of Voice Output	Synthesized speech options. Protective case with Bluetooth amp for high quality speech output.
Language System	The PriO includes LAMP Words For Life which combines the power of UNITY with the principles and strategies of Language Acquisition through Motor Planning (LAMP).
Operating System/Memory	16 GB; iOS 9
Battery Run Time/Charger	Up to 10 hours. Charging via power adapter or USB to computer system.
Access Methods	Wi-Fi and Bluetooth capable. Carry handle, speakers, and kick stand included. Switch access available through single and dual switch ports.
Available Colors	Blue, Black (PriO trim colors)

If the end user of this device has a ventricular shunt, please note that some devices such as the iPad 2<sup>nd</sup> and 3<sup>rd</sup> generation versions contain magnets inside them. There have been studies done indicating that magnets may alter ventricular shunt settings if placed in very close contact (within 10 cm) to the programmable valve of the shunt. You can minimize the risk of this happening by keeping toy magnets, electronic devices, and any other devices containing magnets at a safe distance away from the head. When regularly using any of these devices, keep them from coming in close contact (within 10 cm) to the valve. The original version of the iPad does not contain magnets, so this caution does not apply to first generation iPads.